



ULTRA AVIATION SERVICES, INC.

General Operations Manual

APPENDIX: 1

SAFETY & SECURITY BULLETIN (U-SMS-09)

Revision: Original

Date: 7/1/2015

Bulletin No: 17-007

Date: 3/Jun/17

SUBJECT: LOADING OF EMA (Electrical Mobility Apparatus)

APPLICABLE TO: ALL

When assigned to work a ticket counter and a passenger has an EMA (electrical wheel chair or mobility device) it is your responsibility to:

1- Ask that passenger if the battery on the device is WET or Dry

2- If wet follow procedures for WET battery.

A_ Disconnect Battery. Insulate terminal using tape

B- Place battery in special battery box and follow instruction in box. Remember when working TCX the passenger must provide the special box

C- Fill the EMA baggage tags complete. In some airlines such as TCX you will be required to sign this form

3- If the battery is dry (gel, seal or lithium)

A- Disconnect battery. Before this is done ask the passenger how to do it. Remember this is part of their body and whom better than them to tell you how is properly done

B- Isolate terminal using tape. Note: Some batteries terminal are simply isolated by unplugging them

C- Fill the EMA baggage tags complete. In some airlines such as TCX you will be required to sign this form

IF you are assigned to load an EMA on the aircraft

1- Check with your supervisor where is the EMA to be loaded.

2- Insure that the EMA is ready for transport. Battery disconnected and terminal insulated

3- Check all the information written in the EMA tag

4- Sign EMA tag if you working TCX, for other airlines you might have to sign a WHCR form or a NOTOC

5- SECURED down the EMA inside the belly or the ULD

Remember, when in doubt DO NOT ASSUME, ASK your supervisor or your training department

Follow airlines (AO) SOP regarding Electrical Mobility Apparatus (EMA)

SAFETY FIRST

PREPARED BY:

YOUR SAFETY and Training Department